

DSX Application Software Ver. 2.1.5

Release Note

System Configuration	
OS	Microsoft Windows7 Professional 32bit/64bit SP1
	Japanese or English
CPU	Intel Xeon®E3-1230 Quad-Core3.2GHz or more
RAM	4GB or more
HDD	500G or more ※The capacity to save images, etc. by the DSX application software depends on the free space in HDD.
Graphic	Equivalent to NVIDIA Quadro 600 1GB or more.
Optical Drive	DVD drive (The offline software can be installed to CD drive.)
Monitor Resolution	1920 x 1080

1. Restrictions of Specifications

1-1 Restrictions remaining from the previous version (Ver. 2.1.4)

No	Function	TicketNo	Description	Countermeasures
1	Zoom Change	#4722	If the image acquired with the digital zoom magnification (the status when the [D] icon is displayed on the slider bar for specifying the magnification) is displayed in [vertical shrinking] or [horizontal shrinking] in [dual screen], the position to display the live image is shifted to the upper-left.	If comparing images in the dual screen, do not use the image acquired with the digital zoom magnification.
2	Z Movement	#4692	If Z movement is performed to the Z limit position while setting [Reference Z position] same as [Z limit position], only the message indicating that it reached [Reference Z Position] is displayed, but the message indicating it reached [Z Limit Position] is not displayed.	This is only a display issue and there is no functional problem.
3	Automatic Stitching	#4756	After completing the automatic stitching, the [Stop] button is not returned to the [Acquisition] button infrequently. If the [Stop] button is pressed during this status, the [Acquisition] button cannot be pressed.	This problem is solved by following procedures. (1) Close the automatic stitching screen. (2) Display the [Live Stitching] screen. (3) Press the [Acquisition] button.
4	HDR	#4744	If the HDR image is acquired while the live display stops, even though the fixed parameters for HDR detail settings is set to OFF, the automatic adjustment process of parameters is not executed.	Enable the live display before acquiring the image.
5	Display of Image Information	#4597	If the image acquired by using Ver1.1.5 and enabling HDR (Halation Removal) is opened with the application of this version, HDR=OFF is displayed in the image information display area.	This is only a display issue, and the opened image contains the information of HDR (Halation Removal) : ON.
6	GUI related	#4604, #4724	The top menu of the software (area where [Acquisition], [Measurement] or [Report] buttons are placed at the uppermost area of the screen) is not displayed in some cases.	Press the Windows key on the keyboard, and click the DSX button on the task bar displayed.
7		#4620	If [Acquisition], [Measurement] or [Report] screens are changed consecutively, it may take some time to change the screen.	Wait until the screen changes.
8		#4704	After the Z escape value is entered, the dotted line is displayed vertically.	This is only a display issue and there is no functional problem.
9		#4720	If the image in 3D display is opened while the multi-screen is displayed, the image is not displayed correctly in some cases.	Move the mouse cursor on the image display area, and then move it slightly.
10	GUI related	#4859	If [Acquisition] is selected from [Measurement] in Top Menu very quickly, the screen may not be displayed correctly.	Select [Acquisition] from [Measurement] in Top Menu again.
11		#4861	Immediately after the screen for registering users in Multi Preview is closed by clicking [Cancel], the acquisition screen may not be displayed correctly.	Click the [Measurement] tab or the [Acquisition] tab to display the screen correctly.

12		#4574	At re-login, when the message related to the anti-halation unit is displayed, the drawing of the background is not completed.	The message only appears during drawing the background. There is no functional problem.
13	Touch Panel	#4845	When moving the stage by instructing on the live image with the touch panel, if the next instruction is made before completing the last stage movement, the message indicating "The last operation may not have been executed." may appear.	After selecting [Image Correction]-[HDR], cancel the selection and operate the stage.
14	Observation setting	#4893	If the reproduction of the condition to lower the zoom magnification is repeated twice continuously, the control of the entire application will be Disabled.	Change the objective lens.
15	Step Measurement	#4892	If the image is opened in the order of 3D→2D→3D, the menu shows Plane Measurement, but the icon shows Step Measurement.	Nothing special
16	Report	#4886	If the size of the image opened in logo is large, it takes some time to save the report. ※ Several minutes for several MB.	Please wait for a while.
17		#4891	If the image is zoomed during saving the report, the internal error occurs.	Close the error message
18	Traveling acquisition	#4897	If the [Stop] button is clicked when the traveling acquisition of the first point is completed, the acquired image may not be displayed.	Go to the measurement screen. The acquired image is displayed there. Save it if necessary.
19	Video acquisition	#5006	After Video acquisition was finished in time out, Live observation and Video acquisition cannot be started.	Please restart application.

1-2 Restrictions caused by Upgrade

No	Function	TicketNo	Description	Countermeasures
1				