

PRECiV

Release Notes

PRECiV 2.1.1 Release Notes (April 2024) Build 32039

| ID | Title | Description |
|--------|--|--|
| R-3849 | Movie is added to PRECiV Core and PRECiV Capture | Movie is now a standard feature for all PRECiV packages in PRECiV 2.1.1. |
| R-3848 | Fix the problem that "Manually confirm device status" after acquisition cannot be used | When "confirm device status after acquisition" is selected, the window asking for device status is empty for manual components. |
| R-3847 | Reduce Lens Collision hazard with warning display when exchanging lens type | <p>The current Z-coordinates defined for the acquisition processes, positions in the position list and the defined Z-limit might not be suitable for the inserted lens attachment type.</p> <p>The objective may collide with the sample when the motorized zoom head is moved.</p> <p>A message is added in PRECiV 2.1.1 when the customer is changing from UIS2 -> XLOB -> SXLOB lens attachment type in PRECiV DSX.</p> |
| R-3846 | Fix the problem that the acquired preview images are too large in Best Image. | When acquiring Best Image with PRECiV DSX, the acquired preview images are too large for the thumbnail preview for a few seconds. |
| R-3845 | Zoom tables for SZ61 and SZ51 are updated | Predefined zoom magnification of SZ61 and SZ51 are different from them showing on the zoom drive. |
| R-3844 | Fix the problem that Instant Panorama and EFI doesn't work with reflected light (fine) algorithm | When using Instant Panorama and EFI simultaneously, the EFI information is not compiled properly, using the Reflected Light (fine) algorithm. |
| R-3842 | Fix the problem that after Shadow Contrast image is acquired, the MIX image is monochrome | After acquisition of Best Image with Shadow Contrast mode, the MIX image is not reverted to color mode. |
| R-3841 | Horizontal and Vertical measurements are added to PRECiV Core and PRECiV Capture | When customers are migrating from Stream Basic or Essentials they are missing the Vertical and Horizontal interactive measurements. It is now added to PRECiV 2.1.1 standard Capture and Core packages. |

| ID | Title | Description |
|--------|--|--|
| R-3840 | Improved resolution with DP75 | <p>DP75 with 4K resolution</p> <p>DP75 camera supports a 4K resolution (3840x2160) with PRECiV 2.1.1.</p> <p>In the context of DP75 integration, the limit of EFI has been increased to ~12 MP (4096x3000 pixels).</p> |
| R-3839 | PRECiV DSX: Best Image settings can be adapted to the Zoom head type automatically | <ul style="list-style-type: none"> • Best image settings specific for types of zoom head are available • Best image settings specific for observation modes are available |
| R-3838 | PRECiV DSX: DSX10-SZH is now compatible with Best Image function | The standard zoom head (DSX10-SZH) is now supported with the best image function. |
| R-3837 | Improved Deep Learning performance using ONNX runtime | <p>The Neural Network processing speed is increased when the computer is not fitted with a dedicated GPU.</p> <p>The processing can now be conducted on the CPU directly.</p> |
| R-3836 | Live defect recognition using Neural Network | The Live AI mode is added to PRECiV Observation tab. |
| R-3835 | Adjustments in EFI algorithm calculation | <p>The parameters for rendering of 3D images acquired with EFI are changed.</p> <p>The new algorithm is now taking into consideration all optical components (objective lens and optical zoom) involved.</p> |
| R-3834 | PRECiV DSX: Repeatability assessment of Z calibration | EVIDENT Service personnel can perform Height verification using standard step master sample. |
| R-3833 | PRECiV DSX: Support of Maker Calibration | EVIDENT Service personnel can use the DSX Maintenance Software for performing DSX1000 calibration. |